



EC Declaration of Product Conformity

Vitality Latex Scented (GD001*)

Non-sterile, powder free, latex single use glove

The manufacturer established in the Community declares that the new medical device and PPE described hereafter, Vitality Latex Scented, is in conformity with the provisions of Council Directive 93/42/EEC as amended, national standards EN455-1:2000, EN455-2:2015, EN455-3:2015 and EN455-4:2009 and is Class I self-certified and

is in conformity with the provisions of Council Directive 89/686/EEC as amended and where such is the case, with national standards transposing harmonized standards EN374-1:2003, EN374-2:2003, EN374-3:2003 and EN420:2003+A1:2009

is identical to the PPE which is the subject of EC certificate of conformity no. 7061 issued by SATRA Technology Centre, Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, UK (Notified Body No. 0321)

conforms to the provisions of Council Directive 89/686/EEC Article 11A under the supervision of the Notified Body, SATRA Technology Centre, (Notified Body No. 0321)

Unigloves are a market-leading manufacturer of premium quality, specialist, single use gloves.

We make and supply over 2 billion gloves annually to medical professionals, as well as specialists in the food, automotive, janitorial and care sectors.

 **Unigloves**
Lakeside Park
Neptune Close
Rochester
Kent
ME2 4LT
 01634 726 516
 unigloves.co.uk

Product Information

Glove Care: Store between 5°C and 25°C in a cool dry environment away from direct sunlight.

Food Contact: Vitality Latex Scented gloves have been tested in accordance with EN1186. They are suitable for use with all food categories, in situations of short term repeat contact.

User Information

Cleaning: Gloves cannot be cleaned and reused. Before use, check gloves for physical damage. Do not use damaged gloves.

Fitting and sizing: Only wear products of a suitable size. Products which are either too loose or tight will restrict movement and will not provide the optimum level of protection.

Allergy advice: This product contains natural rubber latex and may contain low levels of residual chemical accelerators, which may cause allergic reactions.

Available sizes

Extra Small	GD0011
Small	GD0012
Medium	GD0013
Large	GD0014
Extra Large	GD0015

Chris Wahlers

Director

Unigloves (UK) Limited



Material safety data sheet

Revision 3 | 24 October 2017

Product code	GD001* The last digit of the order code is dependent on the glove size e.g. 1 = extra small, 5= extra large. N.B. Not all gloves are available in all sizes.	
1. Identification of the substance	Natural rubber latex	
2. Intended application	The intended application of this product are only for the purposes of examination, medical treatments and handling contaminated medical devices to prevent cross contamination between users and the patients. This product is a non-sterile single use natural rubber latex glove and must not be used for sterile applications.	
3. Composition / ingredients	Natural Rubber Latex, ZDBC, ZDEC, ZnO, ZMBT, Sulphur, TiO ₂ , Anti-Oxidant, Wax, Colour Pigment, Vitamin E Acetate, Lanolin Citrus/Peppermint Flavour.	
4. Identification of hazards	This product contains natural rubber latex and accelerators and may cause allergic reactions.	
5. First-aid measures	In case of rashes or any other discomfort appear, discontinue use immediately and wash hand with plenty of water. Consult physician.	
6. Fire-fighting measures	Extinguishing media: water, foam or dry powder. In case of fire wear respiratory protection equipment.	
7. Accidental release measures	No special procedure required. Sweep into appropriate container for disposal. Accidentally released gloves shall not be reused.	
8. Handling and storage	Store in cool and dry area. Protect from sources of heat, direct sun light and UV light.	
9. Exposure control and personal protection	Discontinue use if rashes or other sign of discomfort occur. Thin cotton gloves may be used as liner to prevent direct skin contact.	
10. Physical and chemical properties	Form:	Elastomer
	Colour:	Green
	Odour:	Inherent
	Surface:	Microroughened
	Solubility in water:	Insoluble
	Soluble in:	None identified
	pH value:	Neutral
11. Stability and reactivity	Stability:	Stable
	Materials to avoid:	Petroleum based solvents
	Conditions to avoid:	Sunlight, UV light and heat
	Hazardous decomposition chemicals:	Carbon and sulphur oxides, aliphatic and aromatic hydrocarbons
12. Toxicology information	Skin irritation:	None identified
	Dermal contact sensitisation:	This product contains natural rubber latex and accelerators and may cause allergic reactions in some individuals
13. Ecological information	Non-biodegradable	
14. Disposal considerations	Landfill; dispose of in accordance with local and national regulations. May be incinerated in a suitable facility using heat recovery	
15. Transport information	Not classified	

Material safety data sheet

Revision 3 | 24 October 2017

16. Regulatory information

Primary regulation:	MDD 93/42/EEC and PPED 89/686/EEC
Classification:	Class I non-sterile medical device / Complex design (Cat III) PPE
Risk phrases:	N/A
Safety phrases:	If it is suspected you may have an allergy, discontinue use of the gloves and seek medical advice.

17. Other information

None

The data shown here is based on current knowledge and experience. The purpose of this safety data sheet is to describe the product in terms of its safety requirements. The data does not signify or imply any warranty with regard to the product properties or performance.

Product specification sheet

Revision 5 | 24 October 2017

1.	Product code	GD001*																																		
2.	Product description	Vitality scented Green Latex (Lano-E)																																		
3.	Material	Latex																																		
4.	Inspection parameters	<table border="1"> <thead> <tr> <th>Criteria</th> <th>Inspection level</th> <th>AQL</th> <th colspan="3">Standard</th> </tr> </thead> <tbody> <tr> <td>Dimensions</td> <td>S-2</td> <td>4</td> <td colspan="3">EN455-2</td> </tr> <tr> <td>Physical properties</td> <td>S-2</td> <td>4</td> <td colspan="3">EN455-2</td> </tr> <tr> <td>Water tight test (1,000ml)</td> <td>G-1</td> <td>1.5</td> <td colspan="3">EN455-1</td> </tr> <tr> <td>Visual defects</td> <td>G-1</td> <td>2.5</td> <td colspan="3">In-house</td> </tr> </tbody> </table>					Criteria	Inspection level	AQL	Standard			Dimensions	S-2	4	EN455-2			Physical properties	S-2	4	EN455-2			Water tight test (1,000ml)	G-1	1.5	EN455-1			Visual defects	G-1	2.5	In-house		
Criteria	Inspection level	AQL	Standard																																	
Dimensions	S-2	4	EN455-2																																	
Physical properties	S-2	4	EN455-2																																	
Water tight test (1,000ml)	G-1	1.5	EN455-1																																	
Visual defects	G-1	2.5	In-house																																	
5.	Glove dimensions	<table border="1"> <thead> <tr> <th rowspan="2">Dimensions</th> <th colspan="5">Size (mm)</th> </tr> <tr> <th>XS</th> <th>S</th> <th>M</th> <th>L</th> <th>XL</th> </tr> </thead> <tbody> <tr> <td>Product code</td> <td>GD0011</td> <td>GD0012</td> <td>GD0013</td> <td>GD0014</td> <td>GD0015</td> </tr> <tr> <td>Palm width</td> <td>75 ± 5</td> <td>80 ± 10</td> <td>95 ± 10</td> <td>110 ± 10</td> <td>≥ 110</td> </tr> <tr> <td>Length</td> <td>≥ 240</td> <td>≥ 240</td> <td>≥ 240</td> <td>≥ 240</td> <td>≥ 240</td> </tr> </tbody> </table>					Dimensions	Size (mm)					XS	S	M	L	XL	Product code	GD0011	GD0012	GD0013	GD0014	GD0015	Palm width	75 ± 5	80 ± 10	95 ± 10	110 ± 10	≥ 110	Length	≥ 240	≥ 240	≥ 240	≥ 240	≥ 240	
Dimensions	Size (mm)																																			
	XS	S	M	L	XL																															
Product code	GD0011	GD0012	GD0013	GD0014	GD0015																															
Palm width	75 ± 5	80 ± 10	95 ± 10	110 ± 10	≥ 110																															
Length	≥ 240	≥ 240	≥ 240	≥ 240	≥ 240																															
6.	Glove thickness	<table border="1"> <thead> <tr> <th>Point of thickness measurement</th> <th>Nominal single wall (mm)</th> </tr> </thead> <tbody> <tr> <td>Cuff (at beaded end)</td> <td>0.09</td> </tr> <tr> <td>Palm (at centre of palm)</td> <td>0.12</td> </tr> <tr> <td>Fingertip (at fingertip)</td> <td>0.13</td> </tr> </tbody> </table>					Point of thickness measurement	Nominal single wall (mm)	Cuff (at beaded end)	0.09	Palm (at centre of palm)	0.12	Fingertip (at fingertip)	0.13																						
Point of thickness measurement	Nominal single wall (mm)																																			
Cuff (at beaded end)	0.09																																			
Palm (at centre of palm)	0.12																																			
Fingertip (at fingertip)	0.13																																			
7.	Physical properties	<table border="1"> <thead> <tr> <th>Parameters</th> <th>Before aging</th> <th>After aging</th> </tr> </thead> <tbody> <tr> <td>Median force at break (N)</td> <td>Min 6.0</td> <td>Min 6.0</td> </tr> <tr> <td>Tensile strength (MPa)</td> <td>Min 18.0</td> <td>Min 14.0</td> </tr> <tr> <td>Ultimate elongation (%)</td> <td>Min 650</td> <td>Min 500</td> </tr> </tbody> </table>					Parameters	Before aging	After aging	Median force at break (N)	Min 6.0	Min 6.0	Tensile strength (MPa)	Min 18.0	Min 14.0	Ultimate elongation (%)	Min 650	Min 500																		
Parameters	Before aging	After aging																																		
Median force at break (N)	Min 6.0	Min 6.0																																		
Tensile strength (MPa)	Min 18.0	Min 14.0																																		
Ultimate elongation (%)	Min 650	Min 500																																		
8.	Proteins	< 50 µg/g extractable latex proteins																																		
9.	Powder residue	Maximum 2.0 mg / glove																																		
10.	Cuff type	Beaded																																		
11.	Design	Ambidextrous																																		
12.	Surface	Microroughened																																		
13.	Colour	Green																																		
14.	Quality assurance	Manufactured under ISO 9001:2008 and 13485:2012																																		

The above data represents the product norm, which may be changed without notice and is intended as a guide only.

Warehouse data sheet

Revision 3 | 25 October 2017

Product code

GD001*

Individual Boxes

Product Code	Size	Barcode (EAN-13)	Weight (g.)	Quantity (gloves)	Dimensions (mm)
GD0011	Extra Small	5060264112012	575.00	100	240 (L) x 120 (W) x 70 (H)
GD0012	Small	5060264112029	625.00	100	240 (L) x 120 (W) x 70 (H)
GD0013	Medium	5060264112036	675.00	100	240 (L) x 120 (W) x 70 (H)
GD0014	Large	5060264112043	725.00	100	240 (L) x 120 (W) x 70 (H)
GD0015	Extra Large	5060264112050	775.00	100	240 (L) x 120 (W) x 70 (H)

Cartons

Product Code	Size	Barcode (EAN-13)	Weight (kg)	Quantity (gloves)	Dimensions (mm)
GD0011	Extra Small	5060264112067	6.15	10 x 100	365 (L) x 250 (W) x 250 (H)
GD0012	Small	5060264112074	6.65	10 x 100	365 (L) x 250 (W) x 250 (H)
GD0013	Medium	5060264112081	7.15	10 x 100	365 (L) x 250 (W) x 250 (H)
GD0014	Large	5060264112098	7.65	10 x 100	365 (L) x 250 (W) x 250 (H)
GD0015	Extra Large	5060264112104	8.15	10 x 100	365 (L) x 250 (W) x 250 (H)

Pallets

Product Code	Size	Volume (m ³)	Weight (kg)	Quantity (cartons)	Dimensions (mm)
GD0011	Extra Small	1.92	536.60	84	1,000 (L) x 1,200 (W) x 1,900 (H)
GD0012	Small	1.92	578.60	84	1,000 (L) x 1,200 (W) x 1,900 (H)
GD0013	Medium	1.92	620.60	84	1,000 (L) x 1,200 (W) x 1,900 (H)
GD0014	Large	1.92	662.60	84	1,000 (L) x 1,200 (W) x 1,900 (H)
GD0015	Extra Large	1.92	704.60	84	1,000 (L) x 1,200 (W) x 1,900 (H)

Containers

Product Code	Size	Container Quantity (cartons)		
		20ft	40ft	40ft High Cube
GD0011	Extra Small	1350	2763	3050
GD0012	Small	1350	2763	3050
GD0013	Medium	1350	2763	3050
GD0014	Large	1350	2763	3050
GD0015	Extra Large	1350	2763	3050